

Track mounting bracket

|<u>/</u>//////

Sensor Specifications

| Sensor Range: | 0.15 to 8 meters (±0.001 m accuracy) |
|---------------------------|---|
| Minimum Distance: | 15 cm - false readings when closer |
| Range Switch Settings: | |
| Cart Setting 🖵 📜 : | For cart or short-range activities |
| People Setting <u>犬</u> : | For people or long-range activities |
| Target Indicator LED OFF: | No target detected. Realign target and try again or use a flat, reflective board (648-07373) as a target to improve detection. |
| Surface Reflections: | Tilt the sensor up 5–10 degrees to avoid reflections from a table surface or the sensor housing. |

Motion Quick Start

The PS-2103 Motion Sensor uses ultrasound to measure position.

Additional Equipment Needed

- PASPORTTM USB Link or Xplorer
- EZscreen or DataStudio[®] software (version 1.5 or later)
- Equipment to mount sensor as needed—ring stand (SE-9451) or Cart Adapter Accessory (ME-6743)

Equipment Setup

- 1. Connect the sensor plug to a PASPORT USB Link.
- Connect the PASPORT USB Link to a USB port on your computer or a USB hub.
 The software launches when it detects a PASPORT sensor. Select a choice from
- The software launches when it detects a PASPORT sensor. Select a choice fro the PASPORTAL window.



800-772-8700 • 916-786-3800 • techsupp@pasco.com • www.pasco.com

012-07258D

Click the Start Toggle Data Open Print Exit to Quit Button to Save DataStudio EZscreen **Record Data** Runs Data Data Graph D Example Position (m) Match 1.9 1.8 1.7 1.6 1.5 Graph 0 Sample Bar Chart Data Position vs. Time Graph 0.5 1.0 1.5 2.0 2.5 3.0 3.5 4.0 4.5 5.0 5.5 6.0 6.5 7.0 7.5 8.0 8.5 9.0 9.5 10.0 Your Score Match 91.4 Score Distance from sensor Match Mode button Match Graph Course buttons

EZscreen Specifications

| EZscreen Range: | 0.15 to 2 meters |
|-----------------------|--|
| Recording Time: | Up to 120 seconds for normal graph use Only 10 seconds for Match Graph activity |
| Scale-to-Fit: | Double-click the Graph to scale the data. |
| Information Tool: | Drag cursor over graph to display X,Y coordinate and slope at a point. |
| Export to DataStudio: | Click the Exit to DataStudio button. |

Motion EZscreen

EZscreen Activity - Match Graph

- 1. Select a Match Graph Course button.
- Set the Range Switch on top of the sensor to the People setting.
 Stand in front of the Motion Sensor and aim the sensor at your midsection (Mount sensor on a support rod if necessary or hold a reflector board for better results). Make sure that you
- can move at least 2 meters away from the Motion Sensor.4. Position the computer monitor so you can see the screen while you move away from the sensor.
- 5. Click the **Start** button.
- 6. Use the brief countdown to align your initial position with the beginning of the Match Graph course. When data recording begins, move your body so the plot of your motion matches the gray Position vs. Time plot shown on the screen.
- 7. Your Score shows how close you followed the Match Graph Course.
- 8. Select another **Match Graph Course** button and repeat, or click the **Match Mode** button to exit match mode.